

 **Silhouette**



Hands-free mobile
telephone solution for
hearing aid users

 **Phonic Ear**[®]

Connected to life

SILHOUETTE - now you can hear your conversations on the mobile telephone clearly and distinctly and at the same time speak hands-free




Imagine...

You have an important meeting and wish to call to say that you are late because of the traffic, but it is difficult. You don't have a hands-free set, and there is a lot of noise.

If you have experienced problems with hearing your conversations on the mobile telephone or if you don't have a hands-free set which can be used with your hearing aid(s), then SILHOUETTE is the solution for you.

Hear conversations clearly and distinctly

SILHOUETTE is a discreet, hands-free behind the ear (BTE) inductor providing you with clear and distinct sound from your mobile telephone. SILHOUETTE is connected to your mobile phone, the ear hook is placed behind your ear next to the hearing aid, and the hearing aid is adjusted to T-position (telecoil). You are now ready to make a call and have a clear relaxed conversation.



Your conversation is transferred via SILHOUETTE directly to your hearing aid, and disturbing noise is blocked.

SILHOUETTE allows you to place your mobile telephone away from your hearing aid during a conversation. This eliminates feedback (high pitched noise), which is often experienced by hearing aid wearers when using a mobile telephone.



Choose between two versions

SILHOUETTE is available in two versions. You can choose between SILHOUETTE with one ear hook or with two ear hooks for both ears.



